

BFX® / CFX® THR Implant Selection Chart

Micro/Nano Hip Instrument Set



	Stem		Cup	
Nano Breeds	CFX		CFX 6mm I.D.	
	#1 +0	One piece stem/head Maximum patient weight: 20 lbs (9.0 kgs)	10mm	
	#1 +2			
	#1 +4			
	Stem	Head	Cup	
Micro Breeds	CFX		CFX 8mm I.D.	
	#2	Maximum patient weight: 40 lbs (18.2 kgs)	8mm	12mm
	#3			14mm
				16mm

Universal Hip Instrument Set



	Stem		Head	Cup		
Small Breeds	CFX	BFX	12mm	CFX 12mm I.D.	BFX	
	#4	#4		12mm +0 12mm +3	18mm	20mm 22mm
	Maximum patient weight: 60 lbs (27.3 kgs)		Maximum patient weight: 30 lbs (13.6 kgs)			
Medium & Large Breeds	CFX	BFX	13mm	CFX 13mm I.D.	BFX	
	#4/#5	#5		13mm +1 13mm +3 13mm +5	19mm	20mm
	#4 stem/#5 neck Maximum patient weight: 80 lbs (36.4 kgs)				<i>additional 13mm trial heads and cup impactors required</i>	
					14mm I.D.	
	#5	#6	14mm	14mm +1 14mm +3 14mm +5	20mm	22mm
	#6	#7				
	#7	#8	17mm	17mm +0 17mm +3	17mm I.D.	24mm
	#8	#9		17mm +6 17mm +9		23mm 25mm 27mm
	#9	#10		17mm +11 17mm +13	29mm 31mm	32mm 34mm
	#10	#11		<i>special order</i>		
		#12				
			22mm	22mm +0 22mm +3 22mm +6	22mm I.D.	30mm 32mm 34mm
<i>additional 22mm trial heads and cup impactors required</i>						

Revision / Hybrid BFX Cup

For use with competitors' stem and head



additional 16mm or 19mm cup impactors required

Cup
16mm I.D.
24mm
26mm
28mm
19mm I.D.
30mm
32mm
34mm